Simplicity™


Hands-Free Pumping Bra Kit

Instructions for
Use and Tips

Pump hands-free
and boost milk
supply

See our instructional videos on YouTube at
LactaMed Simplicity™ Bra

View an instructional video at
ittybittyurl.com/1ndb

Neck spring clip at
top of pumping strap
is used to adjust the
strap for an individual fit.

Pumping strap slips
over the head.


Cushioned neck strap
provides extra support
around the neck.

Elastic waist strap with alligator
clip secures the Simplicity™ bra
around the torso.


Spring clips on ends of
pumping strap hold breast
pump shields in place.

Adjusting your Simplicity™ bra for first-time use

Follow these simple steps to adjust the Simplicity™ bra to fit your body the first time you use it. After adjusting it once, you will not need to adjust the fit again.

1. Spread the bra on a flat surface to identify the parts. The pumping strap is designed to be worn around your neck and to support the breast pump shields. The elastic waist strap is designed to circle your torso and secure the bra in place. The clips on each strap are used to adjust the fit of the Simplicity™ bra to each woman's body. Locate the neck spring clip at the top of the pumping strap.
2. Slip the pumping strap over your head so that the neck spring clip rests against the back of your neck. You are now wearing the bra like a necklace.
3. Pressing firmly on the neck spring clip, adjust the length of the pumping strap so that the pumping supports (the ends of the pumping strap) are level with the bottom of your breasts.

4. Wrap the elastic waist strap around your torso and secure it in place using the alligator clip. The waist strap should be near the bottom of your breasts and tight enough to keep the Simplicity™ bra in place.
5. Slip a breast pump shield through one of the pumping support loops at the ends of the pumping strap. Press firmly on the spring clip above the shield and slide it down the pumping strap to secure the breast pump shield securely in place. Repeat on the other side.
6. Position the breast pump shields over your nipples. You may need to slide the pumping strap along the waist strap to adjust the position. You may also need to readjust the length of the pumping strap by using the neck spring clip. Take care to adjust the length and position of the straps carefully to fit your body. Once you achieve the right fit, you never need to adjust the straps again.

Hints:

- Hand Wash – Line Dry.
- A second roll of elastic may be included in your kit. Use this extra strap on top of the breast shields, or as a replacement waist strap.
- If you experience discomfort in the neck area, use cushioning material under the pumping strap or readjust the pumping strap to loosen it. If discomfort persists, discontinue use of the product.

Using Simplicity™ with hands-on pumping (breast massage while pumping with an electric breast pump) can boost your daily milk supply and make pumping more efficient.*

Hands-on Pumping: How To Do It

Most electric breast pumps, even the most expensive ones, are inefficient at emptying the breast of milk. By using the Simplicity™ Hands-Free Pumping Bra kit in combination with hands-on pumping, you can pump more quickly and also increase your daily milk flow.* Follow these simple steps:

1. Use your Simplicity™ bra to hold your breast pump shields in place. Your hands are now free to massage your breasts while pumping.
2. While the electric breast pump is pumping, use your thumbs to massage each breast in long strokes, beginning at the outer edge of the breast and stroking toward the nipple, gently but firmly. Your goal is to encourage the milk stored in the milk ducts throughout the breast to flow toward the nipple.
3. Continue this stroke around the breast, as if you are tracing the spokes of a wheel.
4. As you massage, notice the sprays of milk that occur in the tunnel of the breast shield. These sprays indicate areas of the breast that might not be emptying fully with the electric breast pump.
5. Massage these areas using long strokes.
6. Pump while massaging until you have fully emptied each breast.

*Source: J. Morton et. al., Combining hand techniques with electric pumping increases milk production in mothers of preterm infants, *Journal of Perinology* (2009) 29 757-764.

Item # 70.0001 (Ver. 4.0)

Collecting Colostrum:

Special tips for mothers of newborns

Colostrum is the first kind of breastmilk a mother produces just before giving birth. Colostrum contains important antibodies that protect a newborn against disease. Mothers of infants who are unable to breast-feed immediately after birth can still collect colostrum to feed their babies. Here are some important tips:

- Breast pumps do not efficiently remove colostrum. Colostrum is best hand expressed from the breast using a massage technique. You can also use a combination of breast pump and massage technique to express colostrum and stimulate your milk supply.
- Express colostrum by hand for the first time within one hour of birth.
- Hand express colostrum five times or more per day to feed your baby and stimulate your milk supply.
- Continue to hand express colostrum until your milk comes in, and then begin to use an electric pump in combination with hands-on massage to express and collect breast milk.

Call us at 866-322-2471,
or email us at customerservice@lactamed.com
Instructions in additional languages available at
<http://www.lactamed.com/support-materials/>

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